

T5/T5^{XR}

ADAPTIVE EQUIPMENT

Transform your life.



NuStep worked closely with physical therapists, healthcare professionals, and home users to develop adaptive equipment and accessories that make the T5 and T5^{XR} safer, easier to use, and more effective for a wide range of users and conditions.

NUSTEP'S ADAPTIVE EQUIPMENT SUPPORTS USERS WITH:

- Stroke
- Upper/Lower body muscular deficiencies
- Spinal cord injury
- Traumatic brain injury
- Multiple sclerosis
- Cerebral palsy
- Parkinson's disease
- Obesity
- Many other conditions



Leg Stabilizer

Offers advanced support for users with lower body deficiencies that require bilateral or unilateral leg and foot alignment.

- Stabilizes legs from abducting or adducting
- Adjustable to keep legs in desired alignment
- Attaches on thigh to minimize knee torque
- Easy attachment and removal
- Lightweight aluminum construction
- Washable neoprene attachment strap, no latex
- Sanitary and cleanable
- Made in the USA



WellGrip™

Helps users with upper body deficiencies grip the NuStep handles and maintain a neutral wrist position during their workouts, ensuring good upper body exercise.

- Maintains stable, correct hand position during workout
- Supports wrist and prevents wrist flexion
- Allows full, natural range of upper body motion
- Several docking points for up to 60° of hand rotation
- Straps securely to the T5/T5^{XR} hand grips
- Easy attachment and removal
- Breathable cushioning material



Transforming Lives®

Order Today

Monday-Friday 8am-5pm EST
1.800.322.2209

www.nustep.com



Foot Secure System

The foot secure system combines the ease of ski boot style adjustment with the effectiveness of two foot straps to secure the feet.

- Adjustable to fit any foot
- Easy ratchet adjustment
- Dual straps for stability
- Configurable heel support
- Can be stored under pedals
- Standard on T5^{XR}



Storage Bag

Store NuStep equipment, clipboards or other items in the storage bag conveniently located on the back of the NuStep seat.

- Durable polyester
- 900+ cu. in. storage space
- Loops for Leg Stabilizers
- Sanitary and cleanable
- Also sold as part of a bundle package



T5^{XR} Oversized Seat

To enhance comfort and support, NuStep offers a 22" wide seat with corresponding arm rests and contact heart rate handles.

- Option only for T5^{XR} model
- Pre-installation required
- 4" total additional width
- Comfort for larger users
- Compatible with accessories



NuStep Transporter

The innovative Transporter helps users move any NuStep recumbent cross trainer into almost any desired location.

- Durable steel construction
- Adjustable handle height
- Easy grip handles
- 4" urethane wheels



Polar[®] Transmitter Belt

The Polar[®] chest strap transmits a user's heart rate to the display console for continuous monitoring.

- Fastens around chest
- Comfortable to wear
- Wireless transmission
- Non-interference coded
- ECG accurate



Floor Mat

This practical floor mat fits under any NuStep to protect flooring and add stability during use.

- Rubber material
- Highly durable
- Rolls out easily
- Fits under any NuStep

Special Bundled Packages

TOTAL ESSENTIALS

- 2 Leg Stabilizers (Bilateral)
- 1 WellGrip Set
- 1 Foot Secure System*
- 1 Polar[®] Transmitter Belt
- 1 Accessory Bag

STABILITY ESSENTIALS

- 1 Leg Stabilizer (Unilateral)
- 1 WellGrip Set
- 1 Foot Secure System*
- 1 Polar[®] Transmitter Belt
- 1 Accessory Bag

**The Foot Secure System comes standard on the T5^{XR}*

