# INTRODUCING THE RB8 RECUMBENT BIKE



### **OUR NEWEST ADDITION: A PREMIUM RECUMBENT BIKE**



> RB8<sub>PRO</sub>: 500 lbs. User weight capacity with 18" wide seat

> RB8<sub>MAX</sub>: 600 lbs. User weight capacity with 22" wide seat

### HIGH-DEFINITION DUAL TOUCH-SCREEN DISPLAY

#### Featuring Real-Time Biofeedback



> Engaging programs

> Easy-to-track metrics



> Balanced Power: Captures real-time power and symmetry

## **Specs** (Subject to change)

#### **Unit Dimensions:**

> Approximately 70" L x 32" W

#### Unit Weight:

> Approximately 250 lbs.

#### **User Weight:**

- > RB8<sub>PBO</sub> accommodates users up to 500 lbs.
- > RB8<sub>MAX</sub> accommodates users up to 600 lbs.

#### Fits users:

> From 4'6" to 6'5"

#### Power:

> 24V DC power supply included.

Warranty: Limited 10-year frame, 3-year parts and 1-year labor warranty (US).

## **Build your RB8**

Part No. Description

22101 RB8<sub>PBO</sub> with standard 18" swivel seat with 360° access

22103 RB8<sub>MAX</sub> with wide 22"

swivel seat with 360° access

**Optional Add-ons** 

Part Number Description

52421 Lap Belt – Accessory

(must be selected at time of order)

59095 Torso Belt – Accessory

(must be selected at time of order)

Accessories

**Part Number** Description

05300 NuStep Transporter - Accessory Kit

41815 NuStep Floor Mat

## ORDER NOW!

Call 800.322.2209 or email sales@nustep.com

