

INTRODUCING THE **UE8_{PRO}** AND THE **UE8_{MAX}**



FIRST UPPER BODY ERGOMETER WITH SWIVEL DRIVE TRAIN



Dual HD Touch Screen Monitor
Intuitive software: Engaging programs & easy to track metrics

12-degree seatback recline

Bidirectional operation

Compatible with WellGrip® accessory

Sit, stand or direct wheelchair access without removal or attachment of the seat

Low starting force with 15 levels of resistance

UE8_{MAX}

Maximum Comfort and Adjustability

> The UE8_{MAX} is the only upper body ergometer that accommodates users up to 600 lbs. with a seat that swivels and reclines.

Integrated Sensor Technology Captures Real-Time Power and Symmetry





Now Taking Orders

Specs

Machine Footprint:

- > Straight Only 74" x 32"
- > Right or Left Use Configuration: 74" x 54"
- > Both Right and Left Configuration: 74" x 76"

Unit Weight:

- > UE8_{PRO}: 310 lbs.
- > UE8_{MAX}: 320 lbs.

User Weight:

- > UE8_{PRO} accommodates users up to 500 lbs.
- > UE8_{MAX} accommodates users up to 600 lbs.

Fits users:

From 4'6" to 6'7" in height.

Power:

Medical grade 24V DC power supply included. Specs subject to modification.

Warranty The UE8_{PRO} and UE8_{MAX} are covered under a limited 10-year frame, 3-year parts and 1-year labor warranty (US).

Build your UE8

Base Model Part No.

20001

Description

UE8_{PRO} with standard 18" swivel seat with 360° access

20003

UE8_{MAX} with wide 22" swivel seat with 360° access

Optional Add-ons*

Part Number

52421

Description

Lap Belt - Accessory
(must be selected at time of order)

59095

Torso Belt - Accessory
(must be selected at time of order)

Adaptive Equipment

Part Number

50685

Description

WellGrip® Set - Accessory Kit

50680

WellGrip® – Left Only

50681

WellGrip® – Right Only

20010

UE8 Adaptive & Accessory Essential Bundle (WellGrip set, hip belt, torso belt, floor mat, transporter, accessory bag)

20011

UE8 Adaptive Essential Bundle (WellGrip set, hip belt, torso belt)

Accessories

Part Number

05300

Description

NuStep Transporter – Accessory Kit

41815

NuStep Floor Mat

50780

Accessory Storage Bag – Accessory Kit

20012

UE8 accessory bundle (floor mat, transporter, accessory bag)