INTRODUCING THE UE8_{PRO} AND THE UE8_{MAX}



FIRST UPPER BODY ERGOMETER WITH SWIVEL DRIVE TRAIN



UE8_{MAX}

Maximum Comfort and Adjustability

> The UE8_{MAX} is the only upper body ergometer that accommodates users up to 600 lbs. with a seat that swivels and reclines.

Integrated Sensor Technology Captures Real-Time Power and Symmetry







Now Taking Orders for Shipment Beginning in February 2021

Specs

Machine Footprint:

- > Straight Only 74" x 32"
- > Right or Left Use Configuration: 74" x 54"
- > Both Right and Left Configuration: 74" x 76"

Unit Weight:

Total system below 300 lbs.

User Weight:

- > UE8_{PRO} accommodates users up to 500 lbs.
- > UE8_{MAX} accommodates users up to 600 lbs.

Fits users:

From 4'6" to 6'7" in height.

Power:

Medical grade 24V DC power supply included. Specs subject to modification.

Warranty The UE8_{PRO} and UE8_{MAX} are covered under a limited 10-year frame, 3-year parts and 1-year labor warranty (US).

Build your UE8

Base Model Part No.

20001 UE8_{PBO} with standard 18"

Description

swivel seat with 360° access

20003 UE8_{MAX} with wide 22"

swivel seat with 360° access

Optional Add-ons*

Part Number Description

52421 Lap Belt - Accessory

(must be selected at time of order)

59095 Torso Belt - Accessory

(must be selected at time of order)

Adaptive Equipment

Part Number Description

50685 WellGrip® Set - Accessory Kit

50680 WellGrip® – Left Only 50681 WellGrip® – Right Only

Accessories

Part Number Description

05300 NuStep Transporter – Accessory Kit

41815 NuStep Floor Mat

50780 Accessory Storage Bag – Accessory Kit

Copyright NuStep[®] LLC 2020 • NuStep[®] is a registered trademark of NuStep, LLC, Ann Arbor, Michigan UE8_ss_112020_v1