



Space Requirements by Model

| Model | Footprint | Weight | Minimum Door & Hallway Width | Outlet Required |
|-----------|-------------|----------|------------------------------|-----------------|
| T4 models | 60"L x 27"W | 210 lbs. | 29" | No |
| T5 models | 73"L x 30"W | 298 lbs. | 32" | No |
| T6 models | 73"L x 30"W | 289 lbs. | 32" | Yes |

Going Up or Down a Flight of Stairs

Please note that if the location for your NuStep is located either up or down a flight of stairs, delivery **cannot be completed** if a chairlift or chairlift rails are present.

In the Room

Ensure the flooring is a flat, stable surface. Confirm there is enough space to safely get on and off the machine, and additional space for those accommodating wheelchair access.